

SAMPLE PRINTOUT – JONES SATELLITE/BASE STATION

Jones Medical Instrument Co.  
200 Windsor Dr. Oak Brook, IL

03/07/2017 12:44 pm BTPS @ 25° C  
999999910 Name = \_\_\_\_\_  
Age=58 Height=68 Sex=F Race=Caucasian  
Jones Medical Inst. Co., SA-3, v04.09



Diagnostic Data:

NHANES III	Pred	Test 2	%Pred	Test 1	%Pred	Test 3	%Pred
FVC	3.75	3.39	90	3.22	86	3.18	85
FEV1	3.01	2.64	80	2.45	80	2.26	75
FEV1%	81	78	96	76	93	71	89
FEF25-75	2.68	1.70	63	2.31	86	1.77	66

Interpretative Data:

Lung function normal  
Best: FVC=3.39 FEV1=2.64 FEV1%=78 FEF25-75=1.77

Graphic Data:



This is an example of a typical FVC test report from the Satellite/Base Station. The Satellite can also perform FVL and Post-Bronchodilator testing and display up to 20 different spirometric parameters. Virtually every aspect of the report can be customized to suit the specific needs of the user.