

Spirometry Testing for Inhaled Insulin – Easy As 1, 2, 3

If you're one of the many providers now considering spirometry testing to qualify and monitor your diabetes patients for Afrezza®, selecting the right spirometer is an important factor. The spirometry testing requirements for Afrezza® include both initial qualification and ongoing monitoring. Consequently, selecting a spirometer that is easy, accurate, maintenance-free, and reimbursable should be a priority.

1 Afrezza® Labeling Requirements

Although the Afrezza® labeling only specifies the FEV1 measurement, it's important to know that at a minimum you need both the actual FEV1 measurement and the FEV1 percent-of-predicted value. Here's what you need and why...

FEV1 (i.e., 4.50L) — Must be recorded at baseline, 6 months later and annually thereafter. A decline of $\geq 20\%$ from baseline should elicit consideration of discontinuing Afrezza®.

FEV1 Percent-of-Predicted (i.e., 84%) — Indicates whether the subject has lung disease when determining qualification for Afrezza®.

Jones spirometers include a convenient Spirometry Tracking Chart that streamlines spirometry scheduling and measurement tracking.

2 Qualifying for Spirometry Reimbursement

To qualify for spirometry reimbursement, a spirometry report must be generated and satisfy the Medicare Description for the applicable spirometry CPT code (i.e.94010). Here is what you need...

Graphic Record (i.e., Time-Volume and Flow-Volume printout)

Total & Timed Vital Capacity (i.e., FVC and FEV1)

Expiratory Flow Rate (i.e., FEF 25-75)

Spirometers that don't meet all of these criteria (i.e., spirometers that only measure FEV1 and/or FEV1 percent-of-predicted) do not qualify for Medicare reimbursement. Jones spirometers conveniently calculate and display all of these measurements and graphics from a single exhalation.

3 Choosing the Right Spirometer

Select a spirometer that meets the labeling requirements, reimbursement requirements, and your budget. We provide you with two convenient options:

The Satellite/Base Station Spirometer – A completely dedicated spirometer with docking station and built-in printer.

The Satellite Plus Spirometer – An economical spirometer that allows you to print reports through your PC.

For complete product information and technical specifications, visit www.jonesmedical.com.

Whichever system best suits your needs, you can rest assured that both are simple, accurate, maintenance-free, and fully reimbursable allowing you to easily and confidently qualify and monitor your inhaled insulin patients for years to come.

Only Jones spirometers include a full 3-Year Parts and Labor Warranty, a 30-Day Money Back Guarantee, and are designed, manufactured and supported in the USA. Choosing the right spirometer company is as important as choosing the right spirometer. We've been supporting the spirometry needs of health care providers longer than anyone in the world.

Put our experience and expertise to work for you. Call us today and you can begin providing spirometry testing for your patients tomorrow.

This information is not a substitute for medical advice and does not provide a complete summary of Afrezza® or any inhaled insulin product or address all considerations in making treatment decisions regarding their use. Always refer to the complete product information provided by the manufacturer.

Afrezza® is a registered trademark owned by MannKind Corporation.

Jones Medical Instrument Company is not sponsored by nor has any affiliation with Sanofi US, Inc., MannKind or the Afrezza® trademark.



SPIROMETRY STARTS HERE®

Jones Medical Instrument Company
200 Windsor Drive • Oak Brook, IL 60523 • (800) 323-7336
www.jonesmedical.com



Satellite/Base Station³ Spirometer

- A completely dedicated spirometer with docking station and built-in printer
- Handheld cordless Satellite weighs only 1 pound, eliminating tubes, cables, etc.
- Extremely easy-to-use
- The Base Station recharges the Satellite when docked
- The built-in printer provides a complete printout in 8 seconds
- Fully Reimbursable
- Full 3-year parts and labor warranty
- Maintenance-Free
- Designed, manufactured and supported in the USA



Satellite³ Plus Spirometer

- An economical spirometer that allows you to print reports through your PC
- Handheld cordless Satellite weighs only 1 pound, eliminating tubes, cables, etc.
- Extremely easy-to-use
- Offers seamless printing via Windows compatible PC for 8.5"x11" reports and EMR compatibility
- Fully Reimbursable
- Full 3-year parts and labor warranty
- Maintenance-Free
- Designed, manufactured and supported in the USA

©Jones Medical Instrument Company. All Rights Reserved.