


**NOW AVAILABLE**

# Now get full-page (8.5 x 11) printouts and easy EMR connectivity with the Satellite spirometer!



SMARTPRINT<sup>APC</sup>  
 SMARTPRINT<sup>APC</sup>  
 SMARTPRINT<sup>APC</sup>  
 SMARTPRINT<sup>APC</sup>  
 SMARTPRINT<sup>APC</sup>  
 SMARTPRINT<sup>APC</sup>  
 SMARTPRINT<sup>APC</sup>  
 SMARTPRINT<sup>APC</sup>  
 SMARTPRINT<sup>APC</sup>  
 SMARTPRINT<sup>APC</sup>  
 SMARTPRINT<sup>APC</sup>  
 SMARTPRINT<sup>APC</sup>  
 SMARTPRINT<sup>APC</sup>  
 SMARTPRINT<sup>APC</sup>  
 SMARTPRINT<sup>APC</sup>  
 SMARTPRINT<sup>APC</sup>  
 SMARTPRINT<sup>APC</sup>  
 SMARTPRINT<sup>APC</sup>  
 SMARTPRINT<sup>APC</sup>  
 SMARTPRINT<sup>APC</sup>  
 SMARTPRINT<sup>APC</sup>



**SAMPLE PRINTOUT**

The new SmartPrint<sup>APC</sup> software program allows you to print an entire spirometry report from your Satellite spirometer directly to your PC's printer. APC stands for "All Printer Capable," so all you have to do is load the SmartPrint<sup>APC</sup> software onto your PC, connect the included cable to your PC's USB port and you're ready to go. Simply connect the other end of the cable to the Satellite's serial port and a complete full-page report will automatically print to your printer! It doesn't get any easier. If you have multiple offices, you can utilize your Base Station for printing at one office and SmartPrint<sup>APC</sup> at the other. And, if you're considering transferring your spirometry records to an EMR system now or in the future, SmartPrint<sup>APC</sup> makes it easy.

**Jones—the smartest choice for your spirometry needs.**

- Computer Requirements**
- Windows 7, 8, 10 or 11
  - USB Port

- Satellite Requirements**
- Version 4.09 or greater\*

\*Older versions can be upgraded

- ▼ Only \$245
- ▼ Plug & Play Simplicity
- ▼ EMR Ready
- ▼ No Need to Purchase a Dedicated Printer
- ▼ Multiple Printout Formats
- ▼ View Test Results on Your PC



SPIROMETRY STARTS HERE®

Jones Medical Instrument Company  
 200 Windsor Drive  
 Oak Brook, IL 60523  
 www.jonesmedical.com

**Call us today for details: 1-800-323-7336**